Safety Information

Thank you for choosing Oxa. Read this leaflet and the accompanying documentation carefully before using the Oxa sensor, garment, the Oxa Life App, Oxa, or other related services. Failure to follow the instructions could result in product damage, malfunction, or injury. Keep all safety information and instructions for future reference and provide them to subsequent users of the product.

Visit support.oxalife.com to:

- Access the latest version of all Oxa product documentation.
- Contact Oxa Support and request assistance, instructions, or information.
- Submit any comments or feedback regarding what is good and what can be improved. Oxa values comments from its users.

Reconnect with your body

- Oxa allows you to monitor and influence your nervous system with the aim of improving your wellbeing. Oxa allows you to reconnect with yourself by strengthening the link between your body and mind through conscious breathing.
- Oxa takes you on a journey through breathing exercises. These exercises leverage proprietary biofeedback in real time to help -you learn how to influence your autonomic nervous system.

Important safety information

Oxa is not a medical device. Please note that Oxa is not designed to diagnose, treat, cure, or prevent diseases or medical conditions. The information and guidance in the Oxa journey are provided for informational purposes only and should not be used as an alternative to advice given by physicians or other health professionals. You must always consult a physician if you have any questions regarding a

medical condition or any changes you wish to make to your activity or sleep based on information or guidance from Oxa.

If you do decide to make any changes to your sleep or activity based on information or guidance from Oxa, you do so entirely at your own risk. It is important to always remain sensitive to your body's responses. For example, if you feel pain, fatigue, or discomfort that is unexpected, repeated, or long-term as a result of making changes to your sleep or activity, it is recommended that you speak to a physician before continuing. Information and guidance from Oxa may be interpreted differently if your physiological functions and responses are significantly different than the average physiological functions and responses due to medical conditions or rare natural differences.

The Oxa system uses passive electrodes, which do not pose a health hazard. They are biocompatible and do not cause irritation.

If you experience any redness or skin irritation due to the electrodes being in contact with your skin for an extended period, remove Oxa immediately. Skin rashes may be caused by friction between the garment and the skin. If any symptoms last for longer than two or three days after you stopped using your Oxa product, contact a physician.

If you, or anyone who intends to wear Oxa, have an epileptic condition, consult a physician before you use Oxa. If you experience any of the following symptoms during use – dizziness, altered vision, eye or muscle twitches, loss of awareness, disorientation, involuntary movement, or convulsions – discontinue use IMMEDIATELY and speak to a physician.

Hyperventilation in power breathing is only suitable for people who are generally healthy. If you have any ongoing health conditions — or if you suffer from anxiety — please consult your doctor before practicing hyperventilation. If you are pregnant, don't practice hyperventilation. Long breath holds could harm your baby.

The Oxa sensor complies with Directive 2014/53/EU and

meets the FCC Radio Frequency Radiation Exposure (SAR) standard.

NOTICE: It is splash proof but should not be used while swimming, showering, or bathing. Do not expose the sensor to liquids while it is in use or charging.

⚠ WARNING: Electronic interference. Oxa's technology processes electronic signals which may interfere with implanted cardiac devices, such as pacemakers. If you have an implanted cardiac device, do not use Oxa without guidance from your physician or cardiologist.

⚠ WARNING: Oxa's sensor is a small component, which presents a risk of swallowing and choking to small children. Keep out of the reach of children under the age of 36 months. Seek immediate medical attention if you know or suspect that a child has swallowed the sensor.

⚠ WARNING: When connected to a power source, the device processes electric current and presents a risk of electric shock. Never charge the Oxa sensor while wearing it. Remove Oxa from your body before charging. Disconnect the charger before wearing it.

Using Oxa

To start using Oxa:

1. Pull down the shirt or bra carefully so the electrodes are positioned on each side of the tip of the breastbone. The shirt or bra fits comfortably against the skin but not "too" tight to restrict breathing.
2. Lightly shake or tap the Oxa sensor a few times to activate it.
3. Snap the Oxa sensor into the front of the garment.

To pair your sensor:

1. The sensor should be clipped into your Oxa garment.
2. Bring your Android or iOS device within range and activate Bluetooth®.
3. Open the Oxa Life App.
4. Create an account, or log in.
5. If the sensor LED is not lit or blinking, shake or tap the sensor to wake it from sleep mode. A blue light starts to flash.
6. In the Oxa Life App, tap the Bluetooth icon.
7. Select the Oxa sensor. The Oxa sensor pairs to the Oxa Life App.

If you have paired the Oxa sensor to the Android/iOS device in the device's Bluetooth settings, forget or unpair the device, and start at step 3 to pair your Oxa sensor directly to the Oxa Life App.

Let the app guide you through your first Oxa experience.

If the signal is poor, press and hold it briefly against your skin where the electrodes are positioned. Wet the skin underneath the electrodes with water to improve the quality of the signal. The signal may be poor in cold and dry conditions or if your skin is too dry.

Care instructions

To wash the garment:

1. Remove the sensor before washing the garment.
2. Place the garment in the wash bag (provided) before washing.
3. Use gentle detergent.
4. Wash the garment in a standard washing machine at a temperature of up to 30 °C (90 °F or warm/cold).
5. Lay the shirt flat or hang the garment to air dry. Do not tumble dry the Oxa bra or shirt.

Recommendations

- A quick wash cycle of 30 minutes at 800 rpm or less.
- Wash garment regularly to keep electrodes clean and sticky.
- Use a duster or lint roller between washes to remove any debris on the electrodes.
- Avoid damaging the black surface of the electrodes with your fingernails or accessories such as belts or rings. The garment should still function properly with minor scratches, however major damage to the electrodes can affect the ECG signal.

Important notices to maximize the garment's lifetime

- Do not use bleach or fabric softener.
- Always use a wash bag (provided).
- Do not wash the garment inside out – the electrodes should not be directly exposed so as to avoid wear and tear.
- Do not wring out the garment.
- Do not iron the garment.
- Do not tumble dry the garment. Instead, hang it up to air dry.
- Keep the garment away from fire.

Maintaining the Oxa sensor and app

Oxa releases updates based on user feedback, experience, and to introduce usability improvements. If an update is available, the Oxa Life App will offer you the option to update the sensor firmware. If you decide not to upgrade the sensor firmware, we cannot guarantee proper functionality after six months.

The Oxa battery is rechargeable and lasts for about 25 hours before needing to be recharged. Recharging Oxa takes about 2.5 hours. The LED will be solid red when the battery is charging. The LED will turn green when Oxa is fully charged. Refer to the user manual or support.oxalife.com to understand the color and blinking patterns of the LED status indicator.

NOTICE: Do not use aggressive cleaning substances such as alcohol or bleach to clean the sensor or sensor casing.

Data and privacy

The Oxa Life app complies with the GDPR (General Data Protection Regulation) and the Oxa system was designed with security in mind, keeping the data footprint to a minimum so as to minimize the implications for privacy in turn.

For more information on data and privacy, please refer to our privacy policy, which includes further details on how Nanoleq uses personal data collected through Oxa Life.

oxalife.com/privacy.

Technical information

AVERAGE BATTERY LIFE	25 hours. Battery life varies depending on the type and frequency of use.
LI-POL BATTERY TYPE	Rechargeable, 90 mAh, not replaceable
CHARGER INPUT CONNECTOR	USB-C - USB-A, 5 V=2 A (max.)
CHARGING TIME	150 - 180 minutes
CHARGING CONDITIONS	0 - 45 °C (32 -113 °F)
LED	The LED used in the product is in the visible spectrum and safe for the user and other people nearby.
CONNECTIVITY	Bluetooth Low Energy
FIRMWARE UPDATE	Automatic firmware updates via the Oxa Life app
SECURITY	All measurement data is encrypted directly on the sensor using state-of-the-art cryptographic protocols, such as ECDH and AES, before being sent to your phone. This encryption ensures that your data cannot be read by any unauthorized party.
BATTERY NOMINAL VOLTAGE	3.7 V

BATTERY CHARGING VOLTAGE	4.2 V. Equipment intended to be supplied by an external 5 V DC supply that complies with the IEC 62368-1 standard.
OPERATING TEMPERATURE	0 -45 °C (32 -113 °F)
SENSOR MATERIALS	Housing material PC/ABS, silicone, gold plated pins, and HDPE

Dispose of the sensor (with the internal battery) and garments in accordance with applicable local laws and regulations. Special handling may be required. The sensor is an electronic device and local legislation for recycling electronics should be followed.

Regulatory Information

Declaration of Conformity

CE Hereby, Nanoleq AG declares that this product is in compliance with the Directive 2014/53/EU. The full text of the EU declaration of conformity is available at the following address: www.oxalife.com

UK Declaration of conformity

UK CA Hereby, Nanoleq AG declares that this product is in compliance with the relevant statutory requirements. The full text of the declaration of conformity is available at the following address: www.oxalife.com

Innovation, Science and Economic Development Canada Compliance

This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device. Device contains IC: 11306A-ISPI507

Radio Frequency Radiation Exposure (SAR)

This equipment complies with FCC/IC radiation exposure limits set forth for an uncontrolled environment and meets the FCC radio frequency (RF) Exposure Guidelines in Supplement C to OET65 and RSS-102 of the IC radio frequency (RF) Exposure rules. This equipment has very low levels of RF energy that it is deemed to comply without maximum permissive exposure evaluation (MPE).

FCC Compliance

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions. (1) This device may not cause harmful interference (2) This device must accept any interference received, including interference that may cause undesired operation. Caution: Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. FCC-ID: 2A8RQXA1; Device contains FCC-ID: 2AAQS-ISPI507

Warranty, limitation of liability, intellectual property rights, and terms of use

Refer to oxalife.com for the latest information regarding the warranty, limitation of liability, intellectual property rights, and GTC.

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